

## MUSC Laboratory Supply Request

Facility Name		Requestor Name	
Location Address		Date of Request	
		Date of Need	
Phone			

- \*\*Important Note\*\*** These requested supplies are to be used for specimens sent to MUSC Hospital Lab only. Avoid using these supplies for specimens sent to other reference labs. (i.e. LabCorp, MUSCP-PC(CFC) Lab, Roper, etc). Your cooperation will allow MUSC Hospital Lab to continue providing this service to your location.
- Directions: Indicate quantity of supplies needed. Send the request via email to labcs@musc.edu OR fax to 843-792-4896. **Supply orders are generally filled within 5-7 business days.** If you have any questions or an urgent request, please call Laboratory Client Services at 843-792-0707 for assistance.

Microbiology Collection	Unit of Meas.	Quantity	Blood Collection Tubes	Unit of Meas.	Quantity	Blood Collection Accessories	Unit of Meas.	Quantity
Blood Culture Bottle - Aerobic	Each	___ each	Light Blue (Na Cit) 2 ml	50/Flat	___ each ___ flat	21G Needle w/ preattached holder	100/box	___ each ___ box
Blood Culture Bottle - Anaerobic	Each	___ each	Light Blue (Na Cit) 4 ml	50/Flat	___ each ___ flat	22G Needle w/ preattached holder	100/box	___ each ___ box
Blood Culture Bottle - Pediatric	Each	___ each	Gray (Kox) 2 ml	Each	___ each ___ flat	21G Butterfly w/ assembly	50/box	___ each ___ box
Aptima Kits (Chlamydia, GC, Trich) Urine (yellow)	50/Box	___ each ___ box	Li Hep Green w/ gel, 4 ml	100/Flat	___ each ___ flat	23G Butterfly w/ assembly	50/box	___ each ___ box
Aptima Kits (Chlamydia, GC, Trich) Multitest Swab (orange)	50/Box	___ each ___ box	Gold SST 4 ml	100/Flat	___ each ___ flat	Infant Lancets	100/box	___ each ___ box
E-swab (bacteria culture)	50/Bag	___ each ___ bag	Gold SST 5 ml	100/Flat	___ each ___ flat	Adult Lancets	100/box	___ each ___ box
Universal Viral Transport Media (red cap, pink liquid)	50/Box	___ each ___ box	Red, plastic (clot activator) 5 ml	100/Flat	___ each ___ flat	<b>Pathology/ Cytology</b>	<b>Unit of Meas.</b>	<b>Quantity</b>
<b>Microbiology Media</b>	<b>Unit of Meas.</b>	<b>Quantity</b>	Red, glass (no additive) 10 ml	50/Flat	___ each ___ flat	CytoLyt (Thin Prep) Vials	25/Flat	___ flat
TSA II 5% Sheep's Blood Agar	10/Pack	___ each ___ pack	Lavender EDTA 2 ml	100/Flat	___ each ___ flat	Brush/Spatula	25/Bag	___ bag
Chocolate II Agar	10/Pack	___ each ___ pack	Lavender EDTA 4 ml	100/Flat	___ each ___ flat	Wallach Pappette/ Broom	25/Bag	___ bag
Inhibitory Mold Agar	10/Pack	___ each ___ pack	Pink EDTA 2 ml	100/Flat	___ each ___ flat	Formalin – 40 ml vials	24/box	___ box
Lowenstein-Jensen Slant	Each	___ each	Pink EDTA 6 ml	100/Flat	___ each ___ flat	Michel's Media	Each	___ each
Thioglycolate Broth	Each	___ each	Royal Blue EDTA (purple stripe)	100/Flat	___ each ___ flat	<b>Misc. Supplies</b>	<b>Unit of Meas.</b>	<b>Quantity</b>
Pages Amoeba Saline	Each	___ each	Royal Blue No Additive (red stripe)	100/Flat	___ each ___ flat	Specimen Biohazard Bags	50/pack	___ pack
<b>Urine Collection</b>	<b>Unit of Meas.</b>	<b>Quantity</b>	ACD-A (yellow glass)	50/flat	___ each ___ flat	Plain Pour-Off Tubes (screw cap)	100/pack	___ bag
Sterile Container (urine or stool)	100/bag	___ each ___ bag	Quantiferon Tube (6 ml dk green Li Hep, no gel)	100/Flat	___ each ___ flat	Metal Free Pour-Off Tubes	25/pack	___ bag
Urine Hats	Each	___ each	<b>Blood Collection Microtainers</b>	<b>Unit of Meas.</b>	<b>Quantity</b>	Amber Pour-Off Tubes	25/pack	___ bag
Urinalysis Preservative Tube (Tiger Top)	100/Flat	___ each ___ flat	Amber SST (for Bilirubin)	50/bag	___ bag	<b>Other Supplies (not listed above)</b>		
Urine Culture Preservative Tube (Gray)	100/Flat	___ each ___ flat	Li Hep Green	50/bag	___ bag			
Urine No Additive Tube (Yellow)	100/Flat	___ each ___ flat	Lavender EDTA	50/bag	___ bag			
Urine Transfer Straw	100/Bag	___ each ___ bag	Red (no additive)	50/bag	___ bag			